

Latex Examination Gloves, Powder Free, Non Sterile

Primary Material: Natural Latex Rubber

Latex Content: Latex Rubber

Color: Natural color to white

Size: Extra-Small, Small, Medium, Large and Extra-Large

Design and Feature: Ambidextrous, beaded cuff with finger textured / palm textured finishing, Polymer Coated

PHYSICAL Dimensions

Property	Medical Glove	Standard
Length (mm)	Min 240	ASTM D 3578
Palm Width (mm)	70 ± 10	ASTM D 3578
XS	80 ± 10	
S	95 ± 10	
M	110 ± 10	
L	120 ± 10	
XL		
Thickness – Single Wall (mm)		
Fingers	Min 0.08	ASTM D 3578
Palm	Min 0.08	
Cuff	N/A	

PHYSICAL PROPERTY

Property	Medical Glove	Standard
Elongation at Break (%)		
Before Aging	Min 650	ASTM D 3578
After Aging	Min 500	
Tensile Strength (MPa)		
Before Aging	Min 18	ASTM D 3578
After Aging	Min 14	
Force at Break (N)		
Before Aging	Median 6	EN 455-2
After Aging	Median 6	
Protein Content	≤ 2 mg/glove	ASTM D6124
		EN 455-3
Protein Content	Max 200ug/dm ²	ASTM D5712
		EN 455-3