

## Nitrile Examination Gloves, Powder Free, Non Sterile

**Primary Material:** Acrylonitrile Butadiene Synthetic Rubber Latex

**Latex Content:** Latex Free

**Color:** Blue

**Size:** Extra-Small, Small, Medium, Large and Extra-Large

**Design and Feature:** Ambidextrous, Straight Finger, Finger Textured/Palm Textured, Beaded Cuff, Polymer Coated

### PHYSICAL Dimensions

Property	Medical Glove	Standard
Length (mm)	Min 240	ASTM D 6319
Palm Width (mm)	70 ± 10	ASTM D 6319
XS	80 ± 10	
S	95 ± 10	
M	110 ± 10	
L	120 ± 10	
XL		
Thickness – Single Wall (mm)		
Fingers	Min 0.05	ASTM D 6319
Palm	Min 0.05	
Cuff	N/A	

### PHYSICAL PROPERTY

Property	Medical Glove	Standard
Elongation at Break (%)		
Before Aging	Min 500	ASTM D 6319
After Aging	Min 400	
Tensile Strength (MPa)		
Before Aging	Min 14	ASTM D 6319
After Aging	Min 14	
Force at Break (N)		
Before Aging	Median 6	EN 455-2
After Aging	Median 6	
Residual Powder Content		ASTM D6124
	≤ 2 mg/glove	EN 455-3