

Nitrile Examination Gloves, Powder Free, Non Sterile

Primary Material: Acrylonitrile Butadiene Synthetic Rubber Latex

Latex Content: Latex Free

Color: Blue

Size: Extra-Small, Small, Medium, Large and Extra-Large

Design and Feature: Ambidextrous, Straight Finger, Finger Textured/Palm Textured, Beaded Cuff, Polymer Coated

PHYSICAL Dimensions

Property	Medical Glove	Standard
Length (mm)	Min 240	ASTM D 6319
Palm Width (mm)	70 ± 10	ASTM D 6319
XS	80 ± 10	
S	95 ± 10	
M	110 ± 10	
L	120 ± 10	
XL		
Thickness – Single Wall (mm)		ASTM D 6319
Fingers	Min 0.05	
Palm	Min 0.05	
Cuff	N/A	

PHYSICAL PROPERTY

Property	Medical Glove	Standard
Elongation at Break (%)		
Before Aging	Min 500	ASTM D 6319
After Aging	Min 400	
Tensile Strength (MPa)		
Before Aging	Min 14	ASTM D 6319
After Aging	Min 14	
Force at Break (N)		
Before Aging	Median 6	EN 455-2
After Aging	Median 6	
Residual Powder Content		ASTM D6124
	≤ 2 mg/glove	EN 455-3